

ASTM Method D5769 Quality Control Reference Material without IS (set of 10x2 ml): 14 components

Art. ID LM24-OP-4011645

Unit 10x2 mL

Description

Validity: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2,4,5-Tetramethylbenzene	[95-93-2]	2	%		
	1,2,4-Trimethylbenzene	[95-63-6]	3	%		
	iso-Octane	[540-84-1]	12	%		
	Benzene	[71-43-2]	1	%		
	Ethylbenzene	[100-41-4]	3	%		
	m-Xylene	[108-38-3]	3	%		
	n-Decane	[124-18-5]	12	%		
	n-Dodecane	[112-40-3]	5	%		
	n-Heptane	[142-82-5]	17	%		
	n-Hexane	[110-54-3]	12	%		
	n-Octane	[111-65-9]	17	%		
	Naphthalene	[91-20-3]	1	%		
	o-Xylene	[95-47-6]	3	%		
	Toluene	[108-88-3]	9	%		