

Custom Solution - 6 components in HNO3 5%

Art. ID	LM24-CUS-70994
Unit	1 x 100 mL bottle
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 24 months (originally closed bottle), 12 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium (Na)	[7440-23-5]	1000	mg/L		
	Magnesium (Mg)	[7439-95-4]	500	mg/L		
	Potassium (K)	[7440-09-7]	250	mg/L		
	Calcium (Ca)	[7440-70-2]	1500	mg/L		
	Silicon (Si)	[7440-21-3]	125	mg/L		
	Phosphorus (P)	[7723-14-0]	100	mg/L		