

Native KOVATS-LEE Retention Index Mix 100 µg/mL in Methylene Chloride

Art. ID CIL-ES-5638
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1593 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Toluene	[108-88-3]	100	µg/mL		
	n-Decane	[124-18-5]	100	µg/mL		
	Naphthalene	[91-20-3]	100	µg/mL		
	n-Pentadecane	[629-62-9]	100	µg/mL		
	Phenanthrene	[85-01-8]	100	µg/mL		
	n-Eicosane	[112-95-8]	100	µg/mL		
	n-Tetracosane	[646-31-1]	100	µg/mL		
	Chrysene	[218-01-9]	100	µg/mL		
	Benzo[a]pyrene	[50-32-8]	100	µg/mL		