

## **Tris (4-Isopropylphenyl) Phosphate (unlabelled) 100 µg/mL in Acetonitrile**

Art. ID CIL-ULM-11048-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

### Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter         | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|---------------------------|--------------|---------------------|-------|--------|--------|
|                  | Isopropylphenyl phosphate | [68937-41-7] | 100                 | µg/mL |        |        |