

**EPA Method 545 Internal Standard PDS Mixture (D, 98%) in 1:1 Methanol:Water**

Art. ID CIL-ES-5584  
Unit 2 x 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Uracil (D4,98%)	[24897-55-0]	1	µg/mL		
	L-phenylalanine (D5,98%) )		0,25	µg/mL		