

N,N'-Diphenyl-p-Phenylenediamine (D14, 98%) 100 µg/mL in Nonane (95% Chem. Purity)

Art. ID CIL-DLM-4880-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 95%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N,N'-Diphenyl-p-Phenylene nediamine (D14,98%)	[2687960-60-5]	100	µg/mL		