

EPA Method 545 Internal Native Std PDS Mixture (unlabelled) in 1:1 Methanol:Water

Art. ID CIL-ES-5585
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Uracil	[66-22-8]	1	µg/mL		
	L-phenylalanine	[63-91-2]	0,25	µg/mL		