

Custom Solution - 15 components 10 mg/L each of in Toluene

Art. ID LM24-CUS-71008
Unit 3 x 10 mL ampoule
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	alpha-HCH	[319-84-6]	10	mg/L		
	beta-HCH	[319-85-7]	10	mg/L		
	gamma-HCH (Lindan)	[58-89-9]	10	mg/L		
	delta-HCH	[319-86-8]	10	mg/L		
	2,4'-DDD	[53-19-0]	10	mg/L		
	4,4'-DDD	[72-54-8]	10	mg/L		
	2,4'-DDE	[3424-82-6]	10	mg/L		
	4,4'-DDE	[72-55-9]	10	mg/L		
	2,4'-DDT	[789-02-6]	10	mg/L		
	4,4'-DDT	[50-29-3]	10	mg/L		
	Dieldrin	[60-57-1]	10	mg/L		
	Endosulfan-I (alpha)	[959-98-8]	10	mg/L		
	Endosulfan-II (beta)	[33213-65-9]	10	mg/L		
	Endrin	[72-20-8]	10	mg/L		
	Isodrin	[465-73-6]	10	mg/L		