

VOC - Mix 2 je 200 µg/mL in Methanol

Art. ID LM24-CUS-10179
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dibromochloromethane	[124-48-1]	200	µg/mL		
	Benzene	[71-43-2]	200	µg/mL		