

Custom Solution - 6 components 10mg/L each of in HNO3 2%

Art. ID LM24-CUS-71038
Unit 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months. (6 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	10	mg/L		
	Mercury (Hg)	[7439-97-6]	10	mg/L		
	Arsenic (As)	[7440-38-2]	10	mg/L		
	Selenium (Se)	[7782-49-2]	10	mg/L		
	Platinum (Pt)	[7440-06-4]	10	mg/L		
	Antimony (Sb)	[7440-36-0]	10	mg/L		