

## **Röntgenkontrastmittel - Mix 6 je 10 µg/ml in Methanol/Wasser (1:1)**

Art. ID LM24-CUS-10184  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

Expiration: 18 Months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	lopamidol	[60166-93-0]	10	µg/mL		
	Diatrizoic acid	[117-96-4]	10	µg/mL		
	lohexol	[66108-95-0]	10	µg/mL		
	loxitalamic acid	[28179-44-4]	10	µg/mL		
	lomeprol	[78649-41-9]	10	µg/mL		
	lopromid	[73334-07-3]	10	µg/mL		