

### **Custom Solution - 9 components 1000 µg/mL each of in Methanol**

Art. ID LM24-CUS-71044  
Unit 5 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Pentane	[109-66-0]	1000	µg/mL		
	n-Hexane	[110-54-3]	1000	µg/mL		
	n-Heptane	[142-82-5]	1000	µg/mL		
	n-Octane	[111-65-9]	1000	µg/mL		
	n-Nonane	[111-84-2]	1000	µg/mL		
	n-Decane	[124-18-5]	1000	µg/mL		
	n-Undecane	[1120-21-4]	1000	µg/mL		
	n-Dodecane	[112-40-3]	1000	µg/mL		
	Cyclohexane	[110-82-7]	1000	µg/mL		