

Custom Solution - 3 components in HNO3 2%

Art. ID LM24-CUS-71048
Unit 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	1000	mg/L		
	Antimony (Sb)	[7440-36-0]	1000	mg/L		
	Selenium (Se)	[7782-49-2]	1000	mg/L		