

Custom Solution - 6 components in HNO3 5%

Art. ID LM24-CUS-71050
Unit 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed). (CAPURAN)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Boron (B)	[7440-42-8]	250	µg/mL		
	Tellurium (Te)	[13494-80-9]	250	µg/mL		
	Silver (Ag)	[7440-22-4]	250	µg/mL		
	Silicon (Si)	[7440-21-3]	250	µg/mL		
	Aluminium (Al)	[7429-90-5]	1000	µg/mL		
	Iron (Fe)	[7439-89-6]	1000	µg/mL		