

Gasoline - Benzene Content (Main method: ASTM D 4053 - EN 238 - EN ISO 22854 - ASTM D 6277)

Art. ID ROF-BE-GA-607
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1203 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|-------|--------|--------|
| | Benzene | [71-43-2] | 0,57 | % vol | | |