

DIESEL - Kinematic Viscosity +40°C (Main method: ASTM D 445 - ISO 3104)

Art. ID ROF-VI40-GO-440
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1202 III E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------|------------|---------------------|--------------------|--------|--------|
| | Kinematic Viscosity at 40 °C | | 2680 | mm ² /s | | |