

Custom Solution - 5 components 100 mg/L each of in HNO3 5%

Art. ID LM24-CUS-71062
Unit 125 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed), 12 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Selenium (Se)	[7782-49-2]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		
	Tellurium (Te)	[13494-80-9]	100	mg/L		