

**Custom Solution - 8 components 100 mg/L each of in HNO3 5%, HCl 2%**

Art. ID LM24-CUS-71063  
Unit 125 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

**Description**

Conservation 24 months (originally closed), 12 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Bismuth (Bi)	[7440-69-9]	100	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		
	Platinum (Pt)	[7440-06-4]	100	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Zinc (Zn)	[7440-66-6]	100	mg/L		