

### **Custom Solution - 5 components in HNO3 5%**

Art. ID LM24-CUS-71071  
Unit 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### Description

Conservation 12 months. (6 months after opening)

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Bismuth (Bi)      | [7440-69-9] | 2                   | mg/L |        |        |
|                  | Lithium (Li)      | [7439-93-2] | 2                   | mg/L |        |        |
|                  | Niobium (Nb)      | [7440-03-1] | 2                   | mg/L |        |        |
|                  | Tungsten (W)      | [7440-33-7] | 2                   | mg/L |        |        |
|                  | Titanium (Ti)     | [7440-32-6] | 10                  | mg/L |        |        |