

Custom Solution - ICP Calibration Standard 4 components 1000 mg/L each of in HNO3 2%

Art. ID LM24-CUS-71082
Unit 2 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	1000	mg/L		
	Magnesium (Mg)	[7439-95-4]	1000	mg/L		
	Potassium (K)	[7440-09-7]	1000	mg/L		
	Sodium (Na)	[7440-23-5]	1000	mg/L		