

## **High Si and P iron - chips**

Art. ID                      BAS-BCS-CRM 206/3  
Unit                         100 g (chips)  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	3,17	%		
	Manganese (Mn)	[7439-96-5]	0,72	%		
	Sulfur (S)	[7704-34-9]	0,049	%		
	Phosphorus (P)	[7723-14-0]	1,63	%		
	Nickel (Ni)	[7440-02-0]	0,068	%		
	Chromium (Cr)	[7440-47-3]	0,053	%		
	Vanadium (V)	[7440-62-2]	0,05	%		
	Copper (Cu)	[7440-50-8]	0,1	%		
	Arsenic (As)	[7440-38-2]	0,019	%		