

High Chromium Iron, Block L 40 mm x W 37 mm x T 10 mm

Art. ID	BAS-CRRM 2/1
Unit	block
Deliverydetails	No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	1,92	%		
	Silicon (Si)	[7440-21-3]	1,18	%		
	Manganese (Mn)	[7439-96-5]	1,11	%		
	Phosphorus (P)	[7723-14-0]	0,097	%		
	Silicon (Si)	[7440-21-3]	0,079	%		
	Chromium (Cr)	[7440-47-3]	14,13	%		
	Molybdenum (Mo)	[7439-98-7]	2,44	%		
	Nickel (Ni)	[7440-02-0]	1,61	%		
	Aluminium (Al) total	[7429-90-5]	0,054	%		
	Copper (Cu)	[7440-50-8]	1,59	%		
	Titanium (Ti)	[7440-32-6]	0,07	%		
	Vanadium (V)	[7440-62-2]	0,063	%		