

## **High Chromium Iron, Block L 48 mm x W 42 mm x T 12 mm**

Art. ID	BAS-CRRM 4/2
Unit	block
Deliverydetails	No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,93	%		
	Silicon (Si)	[7440-21-3]	0,45	%		
	Manganese (Mn)	[7439-96-5]	0,58	%		
	Phosphorus (P)	[7723-14-0]	0,049	%		
	Silicon (Si)	[7440-21-3]	0,042	%		
	Chromium (Cr)	[7440-47-3]	21,93	%		
	Molybdenum (Mo)	[7439-98-7]	1,15	%		
	Nickel (Ni)	[7440-02-0]	0,58	%		
	Aluminium (Al) total	[7429-90-5]	<0,005	%		
	Copper (Cu)	[7440-50-8]	0,53	%		
	Titanium (Ti)	[7440-32-6]	0,008	%		
	Vanadium (V)	[7440-62-2]	0,11	%		