

### **Battery Alloy, Block L 50 mm x W 50 mm x T 20 mm**

Art. ID BAS-PBRM L21.01-2  
Unit block  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Copper (Cu)       | [7440-50-8] | 0,0003              | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 0,11                | %    |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | 0,002               | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | <0,0005             | %    |        |        |
|                  | Cadmium (Cd)      | [7440-43-9] | 0,0003              | %    |        |        |
|                  | Arsenic (As)      | [7440-38-2] | <0,0005             | %    |        |        |
|                  | Antimony (Sb)     | [7440-36-0] | 0,0006              | %    |        |        |
|                  | Bismuth (Bi)      | [7440-69-9] | <0,001              | %    |        |        |
|                  | Calcium (Ca)      | [7440-70-2] | 0,051               | %    |        |        |
|                  | Silver (Ag)       | [7440-22-4] | <0,001              | %    |        |        |
|                  | Aluminium (Al)    | [7429-90-5] | 0,012               | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | Bal                 | %    |        |        |