

Battery Alloy, Block L 50 mm x W 50 mm x T 20 mm

Art. ID BAS-PBRM L21.03-2
Unit block
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	<0,0005	%		
	Tin (Sn)	[7440-31-5]	0,38	%		
	Zinc (Zn)	[7440-66-6]	0,002	%		
	Nickel (Ni)	[7440-02-0]	<0,0005	%		
	Cadmium (Cd)	[7440-43-9]	0,003	%		
	Arsenic (As)	[7440-38-2]	<0,0005	%		
	Antimony (Sb)	[7440-36-0]	<0,0005	%		
	Bismuth (Bi)	[7440-69-9]	0,023	%		
	Calcium (Ca)	[7440-70-2]	0,089	%		
	Silver (Ag)	[7440-22-4]	0,01	%		
	Aluminium (Al)	[7429-90-5]	0,011	%		
	Lead (Pb)	[7439-92-1]	Bal	%		