

DL-Isocitronensäure-3Na 500 µg/ml in Methanol

Art. ID LM24-CUS-10188
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Expiration: 18 Months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Trisodium Isocitric Aci d	[1637-73-6]	500	µg/mL		