

## **Standard Cigarette for Ignition Resistance Testing**

Art. ID	NIST-1196a
Unit	2 cartons(400 cigarettes)
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM) is intended for use by laboratories to test the cigarette ignition resistance of soft furnishings and their components, and thermal insulation for resistance to cigarette ignition in accordance with 16 CFR 1632, California Technical Bulletin 117-2013, and 16 CFR 1209. A unit of NIST-1196a consists of two cartons of cigarettes containing 10 packs of 20 cigarettes each. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
on 6.35 mm brass plate	Ignition Strength		95,6 ± 2,0	%	ASTM E2187	
plus 2 layers of filter					, modified	
paper					per NIST	
					TN 1627	