

**ARCONIC - Aluminium base, SUS, Mg, cast, Disc 64 mm x 37 mm**

Art. ID ALCOA-SQ-17  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	~0,005	%		
	Bismuth (Bi)	[7440-69-9]	~0,08	%		
	Chromium (Cr)	[7440-47-3]	~0,25	%		
	Copper (Cu)	[7440-50-8]	~0,35	%		
	Gallium (Ga)	[7440-55-3]	~0,03	%		
	Iron (Fe)	[7439-89-6]	~0,4	%		
	Lead (Pb)	[7439-92-1]	~0,1	%		
	Magnesium (Mg)	[7439-95-4]	~1,6	%		
	Manganese (Mn)	[7439-96-5]	~0,12	%		
	Nickel (Ni)	[7440-02-0]	~0,12	%		
	Silicon (Si)	[7440-21-3]	~0,7	%		
	Tin (Sn)	[7440-31-5]	~0,1	%		
	Titanium (Ti)	[7440-32-6]	~0,08	%		
	Vanadium (V)	[7440-62-2]	~0,03	%		
	Zinc (Zn)	[7440-66-6]	~0,12	%		