

**Custom Solution - 1 component Nitrates as N from Potassium nitrate (N(NO<sub>3</sub>-) from Potassiumnitrate) 100 mg/L in Water**

Art. ID LM24-CUS-71130  
Unit 125 mL bottle  
Deliverydetails No Dangerous Good /not restricted

Description

Conservation 24 months. (12 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Nitrates as N from Potassium nitrate	Nitrates as N (N(NO <sub>3</sub> -))		100	mg/L		