

Yessotoxin (YTX) standard solution in Methanol

Art. ID	NRC-CRM-YTX-c
Unit	0,5 mL
Deliverydetails	De Minimis UN1230 (6.1) II E2

Description

Each ampoule contains approximately 0.5 mL of solution yessotoxin (YTX) at a concentration of 4.3 µmol/L in methanol. The concentration of YTX is suitable for the calibration of LC-MS instruments.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Yessotoxin (YTX)	[112514-54-2]	4,3	µmol/L		