

Polyethylene Resin

Art. ID	NIST-1474b
Unit	60 g
Deliverydetails	No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM) is intended for use in calibration and performance evaluation of instruments used in polymer technology and science for the determination of the melt flow rate using ASTM D1238-13, Test Method for Flow Rates of Thermoplastics by Extrusion Plastometer Standard Test Condition 190/2.16. A unit of NIST-1474b consists of approximately 60 g of white polyethylene pellets in an amber glass bottle. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 190.0 °C ± 0.1 °C	Melt flow rate (FR)		5,01 ± 0,38	g/10 min	ASTM D1238	-13