

**Electrolytic Nickel in Iron, 5 pcs, 4x10x15 mm each (JAPAN)**

Art. ID NMIJ-1006a-1010a  
Unit 5 pcs

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
NMIJ-1006-a	Nickel (Ni)	[7440-02-0]	5,04	%		
NMIJ-1007-a	Nickel (Ni)	[7440-02-0]	10,05	%		
NMIJ-1008-a	Nickel (Ni)	[7440-02-0]	20,02	%		
NMIJ-1009-a	Nickel (Ni)	[7440-02-0]	39,92	%		
NMIJ-1010-a	Nickel (Ni)	[7440-02-0]	60,07	%		