

**Electrolytic Nickel in Iron, 5 pcs, 4x10x15 mm each (JAPAN)**

Art. ID NMIJ-1006a-1010a  
Unit 5 pcs

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| NMIJ-1006-a      | Nickel (Ni)       | [7440-02-0] | 5,04                | %    |        |        |
| NMIJ-1007-a      | Nickel (Ni)       | [7440-02-0] | 10,05               | %    |        |        |
| NMIJ-1008-a      | Nickel (Ni)       | [7440-02-0] | 20,02               | %    |        |        |
| NMIJ-1009-a      | Nickel (Ni)       | [7440-02-0] | 39,92               | %    |        |        |
| NMIJ-1010-a      | Nickel (Ni)       | [7440-02-0] | 60,07               | %    |        |        |