

Aripiprazole-d8 solution, 100 µg/mL in acetonitrile, ampule of 1 mL, certified reference material

Art. ID SAF-A-081-1ML
Unit 1 mL
Deliverydetails De Minimis UN1648 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|----------------------------|----------------|---------------------|-------|--------|--------|
| | Aripiprazole-d8 (Butyl-d8) | [1089115-04-7] | 100 | µg/mL | | |