

**Benzyl piperazine-D7 dihydrochloride solution, 100 µg/mL in methanol (as free base), ampule of 1 mL, certified reference material**

Art. ID                      SAF-B-907-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter                    | CAS number     | Concentration/Value | Unit  | Method | Source |
|------------------|--------------------------------------|----------------|---------------------|-------|--------|--------|
|                  | Benzyl Piperazine-D7 dihydrochloride | [1141738-08-0] | 100                 | µg/mL |        |        |