

Methylphenidate-d9 hydrochloride solution,(Racemic mixture of erythro and threo isomers), ampule of 1 mL, 100 µg/mL in methanol (as free base), certified reference material

Art. ID SAF-M-127-1ML
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Methylphenidate-D9 hydr ochloride (Racemic Mixt ure)	[1219804-02-0]	100	µg/mL		