

Methylone-D3 hydrochloride solution, 100 µg/mL in methanol (as free base), ampule of 1 mL, certified reference material

Art. ID SAF-M-141-1ML
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|----------------------------|----------------|---------------------|-------|--------|--------|
| | Methylone-D3 hydrochloride | [1246820-21-2] | 100 | µg/mL | | |