

**Methylone-D3 hydrochloride solution, 100 µg/mL in methanol (as free base), ampule of 1 mL, certified reference material**

Art. ID                      SAF-M-141-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter          | CAS number     | Concentration/Value | Unit  | Method | Source |
|------------------|----------------------------|----------------|---------------------|-------|--------|--------|
|                  | Methylone-D3 hydrochloride | [1246820-21-2] | 100                 | µg/mL |        |        |