

Ethyl-beta-D-glucuronide-d5, 100 µg/mL, 100 µg/mL in methanol, ampule of 1 mL, certified reference material

Art. ID SAF-E-048-1ML
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------|----------------|---------------------|-------|--------|--------|
| | Ethyl-beta-D-Glucuronid e-D5 | [1305070-98-2] | 100 | µg/mL | | |