

**Ethyl-beta-D-glucuronide-d5, 100 µg/mL, 100 µg/mL in methanol, ampule of 1 mL, certified reference material**

Art. ID                      SAF-E-048-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter               | CAS number     | Concentration/Value | Unit  | Method | Source |
|------------------|---------------------------------|----------------|---------------------|-------|--------|--------|
|                  | Ethyl-beta-D-Glucuronid<br>e-D5 | [1305070-98-2] | 100                 | µg/mL |        |        |