

**Morphine-D3-3-beta-D-glucuronide solution, 100 µg/mL in methanol with 0.05% NaOH, ampule of 1 mL, certified reference material**

Art. ID                      SAF-M-017-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Morphine-3-beta-D-Glucuronide-D3	[136765-44-1]	100	µg/mL		