

**(+/-)-Propoxyphene-D5 solution, 1.0 mg/mL in acetonitrile, ampule of 1 mL, certified reference material**

Art. ID                      SAF-P-904-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1648 II E2

| Text/Information | Analyte/Parameter   | CAS number    | Concentration/Value | Unit  | Method | Source |
|------------------|---------------------|---------------|---------------------|-------|--------|--------|
|                  | (±)-Propoxyphene-D5 | [136765-49-6] | 1                   | mg/mL |        |        |