

**3-Trifluoromethylphenylpiperazine (TFMPP)-D4 hydrochloride solution, 100 µg/mL in methanol (as free base), ampule of 1 mL, certified reference material**

Art. ID                      SAF-T-920-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-Trifluoromethylphenyl piperazine (TFMPP)-D4 h ydrochloride	[15532-75-9]	100	µg/mL		