

**Oxycodone-D3 solution, 100 µg/mL in methanol, ampule of 1 mL, certified reference material**

Art. ID                      SAF-O-005-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number    | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|---------------|---------------------|-------|--------|--------|
|                  | Oxycodone-D3      | [160227-46-3] | 100                 | µg/mL |        |        |