

**Dopamine-D4 hydrochloride solution, 100 µg/mL in methanol with 5% 1 M HCl (as free base), ampule of 1 mL, certified reference material**

Art. ID                      SAF-D-072-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dopamine-d4 Hydrochloride	[203633-19-6]	100	µg/mL		
	de					