

**(1S,2S)-(+)-Pseudoephedrine-d3 hydrochloride solution, 100 µg/mL in methanol (as free base), ampule of 1 mL, certified reference material**

Art. ID SAF-P-056-1ML  
Unit 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	PSEUDOEPHEDRINE-D3   [284665-25-4] Ochloride		100	µg/mL		