

(1S,2R)-(+)-Ephedrine-d3 hydrochloride solution, 100 µg/mL in methanol (as free base), ampule of 1 mL, certified reference material

Art. ID SAF-E-025-1ML
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1S,2R(+)-Ephedrine-D3 h ydrochloride	[285979-73-9]	100	µg/mL		