

(1S,2R)-(+)-Ephedrine-d3 hydrochloride solution, 1.0 mg/mL in methanol (as free base), ampule of 1 mL, certified reference material

Art. ID SAF-E-026-1ML
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1S,2R(+)-Ephedrine-D3 h ydrochloride	[285979-73-9]	1	mg/mL		