

**Phenobarbital-D5 (ring) solution,(deuterium label on ring), 100 µg/mL in methanol, ampule of 1 mL, certified reference material**

Art. ID                      SAF-P-018-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phenobarbital-D5 (Deuterium Label On Side Chain)	[72793-46-5]	100	µg/mL		