

**5alpha-Dihydrotestosterone-D3 (16,16,17-D3) solution, 100 µg/mL in methanol, ampule of 1 mL, certified reference material**

Art. ID                      SAF-D-077-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	5alpha-Dihydrotestoster one-D3 (16,16,17-D3)	[79037-34-6]	100	µg/mL		